# Giardiasi

**IDEXX**  
LABORATORIES

## Indice

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*Giardia duodenalis* (syn. *Giardia intestinalis* and *Giardia lamblia*) is an intestinal flagellate infecting human beings and a wide range of mammals, including dogs. Currently, 8 distinct assemblages or genotypes (A–H) have been identified within this protozoan species. In particular, dogs can be infected by host-specific assemblages C and D or zoonotic assemblages A and B. Infective cysts are present in the feces of infected hosts, and transmission occurs by accidental ingestion of *Giardia* cysts contaminating soil, coat, food, and mostly water.

***Giardia duodenalis* is one of the most common enteric parasites of dogs worldwide.**

Because genetically identical, potentially zoonotic genotypes (predominantly assemblage A) may exist in human beings and dogs sharing the same environment, recent studies have focused on the role of dogs as a source of zoonotic transmission of *G. duodenalis* in human communities. The sampling of dog stools found in the environment could therefore be convenient and serve for screening studies in human settlements. For this purpose, knowledge of appropriate tests for programs dealing with the risk of environmental contamination with *G. duodenalis* is important,

as the effectiveness of these programs depends on the sensitivity (Se) and the specificity (Sp) of the tests used to identify contaminated dog feces in the environment. While Se and Sp are fixed test properties, they can be influenced by the prevalence of *G. duodenalis* and this influence thus affects test accuracy.

In veterinary medicine, the fecal flotation method with zinc sulfate is regarded as the best choice for the microscopic identification of *Giardia* cysts, with the test showing a Sp of **100%**. Due to the intermittent nature of cyst shedding, the Se of this procedure is **70%** when a single sample is evaluated and increases to **95%** if at least 3 fecal specimens are examined within 5 days. Unfortunately, the possibility of multiple sampling from the same subject is often impractical and almost impossible when dog feces found in the environment are being tested.

Other diagnostic techniques have been developed for diagnosis of canine giardiasis in a single test and include identification of cysts by a direct fluorescent antibody test (DFAT), detection of coproantigens by traditional enzyme-linked immunosorbent assay (ELISA), and DNA amplification by polymerase chain reaction (PCR)-based methods. However, such methods are costly and require sophisticated and expensive equipment as well as trained personnel. Therefore, the tests are best suited for well-equipped diagnostic centers. In frontline veterinary clinical laboratories, ELISA-based rapid diagnostic test kits are often used routinely as practical and

time-saving assays. One particular commercial ELISA kits has become available in Italy and many other countries for the detection of *Giardia* antigen in canine and feline feces.

Microscopic examination for diagnosis of canine giardiasis is an “alloyed” gold standard test (i.e., an imperfect reference test) with **100%** Sp but lacking Se. A Bayesian model, which estimates the Se and Sp of 1 or more tests in the absence of a gold standard, can help in identifying the most suitable tests for routine diagnosis and screening programs. This is particularly useful when none of the competing tests are treated as the gold standard, and the diagnostic error rates are estimated for all studied tests. Prior information (such as that generated by previous research or expert opinions) is combined with the observed data to obtain posterior distributions of test parameters.

Knowledge of the true disease or infection status of the animal is therefore not necessary, and instead this unknown information is incorporated into the model as a latent variable. Two studies have been published to date using a Bayesian model to estimate the performance of diagnostic techniques for canine giardiasis, and only 1 of these included evaluation of the commercial kit evaluated in the current study.

A previous study, carried out using a real-time PCR assay, revealed that **30.8%** of dog feces left uncollected on soil in urban areas harbored *G. duodenalis* DNA. The present study represents an extension of the previous survey as the same samples were re-examined using a commercial ELISA with the performance of the ELISA evaluated using a Bayesian model. [...]

It has been reported that *G. duodenalis* may have an impact on public and animal health through environmental contamination with feces of infected hosts. The present Bayesian evaluation provides evidence on the efficiency of the assay for the detection of *Giardia* antigen in dog feces dispersed in the environment as an indicator of field contamination. Based on the results of the current study, the test is more likely to have false-negative results. In a clinical setting, the high Sp minimizes the risk of false positives and unnecessary antimicrobial treatment. On the other hand, a proportion of infected dogs may go undetected due to false-negative results, though to a limited extent. From a public health point of view, infected dogs may serve as reservoirs and be potential sources of infection for human beings and other animals, mostly when they remain undetected and thus untreated. It has been suggested that the Se of the assay can be increased by multiple samplings. This can be done in dogs

but is unrealistic for studies on dog feces found in the environment. If used on canine feces contaminating the environment, as well as on symptomatic dogs with diarrhea, the assay could be performed as a first step for primary screening while another test could be carried out supplementary to the assay, to improve the likelihood of correctly identifying dog feces with *G. duodenalis* cysts. For instance, additional diagnostic testing with PCR-based methods could provide added sensitivity and would be indicated because molecular characterization of isolates could be performed to determine their assignment to genotypes.

The present results show that the assay, which features rapid processing of fecal samples without a complex technical structure and extensive costly labor, is a specific and fairly sensitive technique for the diagnosis of *Giardia* in clinically affected dogs as previously reported. As well, **the current study has shown that the ELISA can also be an appropriate method** and a practical alternative to other methods for large-scale environmental screening programs **to assess the risk of exposure to *G. duodenalis* cysts of canine source in human settlements.** [...]

As the ELISA and the real-time PCR assay measure different biological parameters associated with *Giardia* in dogs (i.e., the ELISA detects the presence of *Giardia* coproantigen while the real-time PCR assay measures *Giardia* cyst DNA in fecal samples), the 2 tests can be considered as independent of each other. Therefore, the ELISA and the real-time PCR assay can be used in combination at the same time in a single animal to allow improved diagnosis. Parallel and series interpretation can be used to combine the information provided by the individual test results. [...]

The present report confirms the **efficiency of the ELISA** for the diagnosis of canine giardiasis and, in addition, provides evidence on its possible application to dog feces found in the field. The main advantages of the assay are that it is rapid, is relatively inexpensive, and employs user-friendly techniques without special equipment. This suggests that the assay could reasonably be included, as a valid alternative or as an adjunct to other techniques, in screening studies for the identification of areas exposed to risk of environmental contamination with *Giardia* cysts of canine source. Nonetheless, molecular procedures for detection and identification of *Giardia* species and genotypes cannot be replaced when the zoonotic potential of canine feces contaminating the environment has to be accurately assessed in human communities.

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© 2013 IDEXX Laboratories, Inc. Tutti i diritti riservati • Tutti i marchi ®/TM sono di proprietà di IDEXX Laboratories, Inc. o di suoi associati negli Stati Uniti e/o in altri paesi. La politica di tutela della privacy di IDEXX è disponibile sul sito [idexx.it](http://idexx.it)
IDEXX Laboratories Italia S.r.l.